
a D²GROUP company

Faster S.r.l.
Italy, Via Merendi, 22
20010 Cornaredo (MI)
TEL. +39 02 93 991 92
FAX +39 02 93 991 608
www.faster-air.it

ISO-9001 Accredited Company

Laminar Airflow - A Guide
Flux d'Air Laminaire - Un Guide

Theoretic principles
Definitions
Classifications
Applications
and practical uses

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TABLE OF CONTENTS

- ✓ A Guide to Terminology
- ✓ Biological risk classification and safety systems
- ✓ Criteria to classify the Microbiological Safety Cabinets
 - Class I Microbiological Safety Cabinet
 - Class II Microbiological Safety Cabinet
 - Class III Microbiological Safety Cabinet
- ✓ General rules for a wise use of the Microbiological Safety Cabinets
 - Working surface
 - Manual skill
- ✓ Cleaning and maintenance of the Microbiological Safety Cabinet
- ✓ Waste treatment and disposal
- ✓ Use of cytostatic drug safety cabinets while preparing therapies with antineoplastic drugs
- ✓ References

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A GUIDE TO TERMINOLOGY

Laminar airflow

A laminar airflow is an unidirectional flow composed of tiny "wisps" of air which, by moving at the same speed in parallel directions, are able to produce a draught without any turbulence.

HEPA filters

A HEPA filter (High Efficiency Particulate Air) is a filter composed of glass microfibres, welded between them by means of an epoxy resin, pleated like an accordion and set within an aluminium frame.

The way of testing the filter provides for an aerosol leakage of DOP (D-2ethylhexylphthalate), of DOS/DEHS (D-2ethylexylbactate) or of Emery 3004 (PolyAlphOlephin). Its efficiency is assessed on the basis of its capacity to contain the 0.3 µm particles produced and more generally must be higher or equivalent to 99.999% DOP (0.001% penetration) in order to ensure, under standard environmental conditions, the class100 (M3.5) according to the U.S. Federal Standard 209/e (equivalent to the ISO 5 class according to the EN ISO 14644-1) and in order to comply with the recent European regulation EN-12469 concerning the Microbiological Safety Cabinets.

The European regulation for the classification of the EN-1822 filters has recently modified the control system of HEPA filters with the introduction of a new method of testing the efficiency: the MPPS (Most Penetrating Particle Size).

Just to provide an example, a HEPA H14 filter has a 99.999% DOP efficiency that corresponds to 99.995% MPPS.

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Class 100 definition

This term is referred to the air particle contamination. The class related to an environment shows in fact the maximum number of particles accepted in a given volume.

The U.S. Federal Standard 209/e is still the most important standard concerning the particle contamination, even if, since November 29th 2001, it has been cancelled. It considers as "class 100" environment, the one containing at most 100 particles whose size varies from 0.5 to 5 μm in a cubic foot volume, nearly 28 litres, or "M3.5" for the same contamination class which is though referred to the European metric measurement system that is to say at most 3530 particles in a cubic metre (1000 litres) of air. The new ISO standard, which replaced the previous U.S. Federal Standard 209/e, called ISO 14644-1, represents the ISO 5 class which is equivalent to the previous name of class 100 (M3.5), with an amount of 100.000 particles of 0.1 μm , 23.700 particles of 0.2 μm , 10.200 particles of 0.3 μm , 3520 particles of 0.5 μm , 832 particles of 1 μm and 29 particles of 5 μm in a cubic metre of air.

The Class 100 vertical and horizontal laminar flow cabinets

We are talking about cabinets which, thanks to fans and HEPA filters, are able to move inside them, a certain amount of laminar or unidirectional flow air and have the contamination characteristics described by the U.S. Federal Standard 209/e relative to class 100 (M3.5) as well as by the EN ISO 14644-1 for the ISO 5 class. As for the horizontal flow type, the air moves from the back wall of the cabinet, where the HEPA filter is placed, going towards the operator. This kind of flow has the advantage of avoiding a possible contamination caused by the operator's hands as well as of creating wide spaces allowing to introduce huge objects and to easily operate. Nevertheless, the operator's protection is not guaranteed.

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As for the vertical flow cabinets, the air coming out of the HEPA filter goes instead downwards, towards the working surface; though allowing a better containment, if we are not using the microbiological safety model, these very cabinets cannot guarantee the operator's protection.

While using the "class 100 laminar flow" cabinets, the air speed is 0.45 m/s +/- 20% according to the specifications of the EU GMP (Volume 4 Pharmaceutical Legislation) as A degree, and while using the Class II microbiological safety cabinets, the European regulation EN 12469:2000 allows the manufacturer to select speed between 0.25 and 0.50 m/s +/-20%.

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BIOLOGICAL RISK CLASSIFICATION AND SAFETY SYSTEMS

When manipulating biological pathogen material or potentially pathogen substances, the most suitable collective safety mean is a microbiological safety cabinet (Biohazard) helping to prevent aero-diffusion of pathogens within the laboratory. The right choice of the Biohazard cabinet must be done by considering the biological risk level of the material that the operator will be handling inside it. Obviously, given that we are talking about "bio-safety requirements", it is essential to draft an Assessment Report about the Risk and to implement some global protocols including IPD, cabinets, laboratory layout, SOP, training and information, diagnostic monitoring of the staff and whatever else is necessary to cancel or to reduce the risk to a minimum. Different international bodies drafted some biological risk classifications, by taking into account both the dangerousness of the pathogens and the minimum requirements to guarantee protection during their manipulation.

The Dangerous Pathogens Advisory Group - Department of Health and Social Security (DPAG, DHSS-GB), the Centre for Disease Control - Public Health Service (CDC-USA), to which also the World Health Organization (WHO) usually addresses, made a list of all the pathogens on the basis of their dangerousness. The National Cancer Institute (NCI USA) classified the oncogenic viruses and the National Institute of Health (NIH-USA) established some parameters to operate with recombinant DNA.

The DPAG, the CDC, the NIH determined 3 or 4 bio-safety levels including a set of approaches, equipment, infrastructure able to guarantee ever increasing levels of protection, both for the staff and the environment.

Besides the product characteristics (pathogenicity level, ways of transmission, resistance to antibiotics, etc..), other parameters to be considered are represented by the amount of material to be handled and the know-how of the working staff.

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The Directive No. 2000/54/CE divided the biological agents in four groups defined as follows.

Hazard Group 1
(it detects low individual and collective risk pathogens)

"A biological agent unlikely to cause human disease"

it includes biological agents such as bacteria or similar organisms, viruses, parasites and mushrooms not showing too much danger for the laboratory staff, as long as all the asepsis rules about microbiological approaches are observed. Generally, we are talking about micro-organisms that are not likely to cause illnesses in men and animals. (WHO - "Bio-safety Manual").

Hazard Group 2
(it detects moderate individual risk and limited collective risk pathogens)

"A biological agent that can cause human disease and may be a hazard to employees; it is unlikely to spread to the community and there is usually effective prophylaxis or effective treatment available"

This risk level includes the biological agents likely to cause diseases in human beings and animals, but they are not likely to represent a serious risk to those working in a laboratory, to the community and to the environment. An accidental exposure in a laboratory can cause infections, but efficient preventive measures as well as treatments are available (WHO - "Bio-safety manual").

The second biological risk level covers the use of a microbiological safety cabinet.

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Hazard Group 3 **(it detects high individual risk and low collective risk pathogens)**

"A biological agent that can cause severe human disease and presents a serious hazard to employees; it may present a risk to spreading to the community but there is usually effective prophylaxis or effective treatment available"

This group includes those pathogens that usually cause serious diseases in men, but that do not usually spread from an affected individual to a healthy one (WHO - "Bio-safety Manual").

In that case, although the microbiological safety cabinets are an excellent primary collective barrier against the infections spreading within the laboratories, it is wise to provide for a further protection barrier in favour of the people living and working outside the laboratory.

With the aim of protecting the external environment against possible leaks of pathogens, some safety laboratories such as BL3 or P3 are installed.

Hazard Group 4 **(it detects high individual and collective risk pathogens)**

"A biological agent causes severe human disease and is a serious hazard to employees; it is likely to spread to the community and there is usually no effective prophylaxis or treatment available"

Handling pathogens belonging to this last group requires the maximum level of biological safety. This group includes micro-organisms, viruses and substances considered as being very dangerous for those who must handle them and representing a high risk of epidemics.

They are not usually used in everyday clinical and research laboratories and their treatment requires the mastery of the rules related to them.

The most suitable laboratory for manipulating such kind of micro-

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organisms must be able to contain the highest amount of biological substances (type BL4 or P4) together with the use of some methods for the staff's protection, such as integral overall provided with pressurized air mask, Class III microbiological safety cabinet and so on.

This kind of laboratory is made up of a type BL3 or P3 lab around which there is a positive pressure air space compared to the environment and the laboratory.

This way, it is possible to create a "dynamic tightness", necessary to complete the "static tightness", which is obtained by sealing each crack in the room.

The laboratory access must be strictly checked and allowed just to authorised staff; a filtered air antechamber (a locker room), equipped with a sterilizing shower by means of chemicals as well as with airtight doors allows operators to wear all the garments and the devices to be used only inside the safety laboratory.

Inside the P4 laboratory, the bare minimum in terms of working tools must be introduced.

The mere primary barriers are not sufficient to prevent, in case of accidents, the spread of pathogens in the environment, thus representing a risk to the whole community. It is necessary, in some cases, above all when handling agents belonging to groups 3 and 4, to operate inside the secondary barriers, that is to say inside the bio-safety laboratories. Contrary to the clean room, where the objective is to safeguard the product sterility by means of a contamination-controlled environment, when we consider a bio-safety laboratory BL2/BL3/BL4, the aim is to safeguard the operator and the community against the risks of a possible contamination by infectious agents.

We can distinguish three levels of containment: the BL2-basis laboratory, the BL3 bio-safety laboratory and the BL4 maximum safety lab.

The BL2-basis laboratory has been planned and fitted with such an equipment to make it ideal to handle biological agents belonging to group 2. It is not possible to use biological agents of a higher group. In short, it must have fittings and finishes allowing to easily clean and to use primary barriers, including the microbiological safety cabinet which shall be placed in the right position. It must also have in its surroundings

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a system for inactivating and disposing of waste in safety conditions.

The BL3 bio-safety laboratory is essential to operate with biological agents of group 3. It is not possible to use higher level biological agents. They are characterized by the presence of an area or entry antechamber (filter zone or locker room) with airtight doors, while inside the laboratory, the necessary primary barriers will be arranged.

Moreover, walls must comply with the pharmaceutical standard, that means they must be smooth and easy to clean, welded between them, all their corners must be rounded off by means of gougues. The incoming air must be treated and filtered by means of a special HVAC (Heating, Ventilation and Air Conditioning) group aiming at providing the best conditions in terms of temperature and humidity as well as of air cleaning. The air rejected by the laboratory is totally filtered, with the cabinets set on HEPA terminal filters, before being ducted outside the building. All the material going out of the safety laboratory must be reclaimed before its disposal.

It can be devised into two ways: a system with "cascade pressure gradients" of the different rooms or with a "positive pressure in the filter zone and atmospheric pressure along the corridor and in the laboratory room" on the basis of the kind of work and to financial means.

The BL4 highest safety laboratory is a very complicated structure, used just in few specialised centres which are authorised to use biological agents belonging to group 4.

In this BL4 lab, the primary barriers used in the previous kinds of laboratories BL2 and BL3 are not sufficient, that's way some special measures must be arranged; for example the operator's protection is guaranteed by pressurized overalls, all activities are carried out inside Class III biohazard cabinets or glove-boxes in order to avoid the direct contact with the products.

Whatever comes out of the laboratory must be sterilised throughout chemical and physical means, properly approved and validated. Given the complexity and the high risk to implement the BL4 laboratories, it is necessary to obtain the authorization by the competent Healthy Authorities.

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CRITERIA TO CLASSIFY THE MICROBIOLOGICAL SAFETY CABINETS

The classification, the description and the technical requirements that the Microbiological Safety Cabinets must have, were established by product standards issued nationwide:

- ✓ NSF/ANSI 49 - 2002
(National Sanitation Foundation - USA)
- ✓ EN 12469 - 2000
(European Standard)
- ✓ AS 2252 - 2004
(Australian Standard)

The microbiological safety cabinets are divided in three different classes, each one with its modes of working and with different structures, able to guarantee different protection levels.

	CABINETS CLASS I	CABINETS CLASS II	CABINETS CLASS III
Operator's protection.....	good.....	good.....	excellent.....
Environmental protection.....	good.....	excellent.....	excellent.....
Product protection.....	none.....	excellent.....	poor.....

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RUNNING INSTRUCTIONS

- A** Front opening
- B** Front screen
- C** Exhaust HEPA filter
- D** Exhaust plenum


- Ambient air
- Contaminated air
- HEPA filtered air

Class I Microbiological safety cabinets

Definition: Class I type safety cabinets provide minimum operator protection but are not designed to protect the product.

These are open-fronted exhaust cabinets inside of which operator can carry out manipulations of low and moderate risk biological agents. Air from the room flowing into the cabinet minimises the escape of aerosols from the work area and is exhausted through the high efficiency air filter system (HEPA). Class I Cabinets are available in both 'ducted' and 'recirculatory' models.

Application: low and moderate risk biological agents.

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RUNNING INSTRUCTIONS

- A** Front opening
- B** Front screen
- C** Exhaust HEPA filter
- D** Main HEPA filter
- E** Positive pressure plenum
- F** Negative pressure plenum


- Ambient air
- Contaminated air
- HEPA filtered air

Class II Microbiological safety cabinets

Definition: Class II MSC-s not only provide the same level of operator protection afforded by Class I MSC-s but will also protect the samples and products on the worktop. In addition to the even and continuous ingress of air through the open front of these cabinets product protection is achieved by the even and unidirectional downward air-flow - referred to as "laminar flow".

The inward flow of air from the room forms an air curtain which minimises the escape of aerosols from the cabinet - providing operator protection. The filtered downward flow of air of a high degree of "laminarity" protects the samples and the work on the worktop - preventing cross-contamination. About 70 % of the air is internally recirculated in this laminar-flow pattern via the HEPA Filter - and about 30 % of the air is exhausted via a second HEPA (high efficiency) air filter.

Application: all the categories of low/moderate risk biological agents.

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RUNNING INSTRUCTIONS


- A** Front with gloves for manipulation
- B** Front screen
- C** Exhaust HEPA filter
- D** Main HEPA filter
- E** Passing sterilizer with a double door or box to pass preparations

- Ambient air
- Contaminated air
- HEPA filtered air

Class III Microbiological safety cabinets

Definition: Class III type Microbiological Safety Cabinets are totally enclosed "glove-box" type safety cabinets - providing total operator protection and are suitable for work with hazard group 4 pathogens. A transfer hatch allows materials and instrumentation to be introduced into the Class III cabinet without breaking the integrity of the sealed enclosure. Gloves fitted upon glove-ports allow the users to manipulate the samples and products - secure in the knowledge of complete isolation from the pathogens and infective or toxic materials being used.

Application: all the categories of biological agents.

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GENERAL RULES FOR A WISE USE OF THE MICROBIOLOGICAL SAFETY CABINETS

Be sure the cabinet is proper to manipulate the sample to be treated and to the specific applications.

Install this equipment in a room fit to its size and far from heat sources, walls, pillars, passages, communicating areas or heat sources that may jeopardize the correct function of the cabinet (consult the instructions written by the manufacturer in the Instructions Manual). Use this cabinet with the motor fans always on and with the front screen properly set and when there are no audible and/or visual alarms.

Duct the expelled air outwards, every time you are handling radioactive tracers or cytostatic drugs. In this case, set an anti-blowback damper in order to avoid possible air recirculates inside the extraction channel. When it is not possible to carry out what mentioned above and when the exhausted air recirculates onto the working place, it is wiser to install a second HEPA absolute filter on the exhaust channel, which can be individually tested or an activated charcoal filter when manipulating chemicals.

Set all the necessary material inside the cabinet, then turn it on and let it run at least 20 minutes before starting the work, and verify there are no alarm signals. While working, check at regular intervals the tools monitoring the safety conditions and pay attention to any alarm, to be sure the cabinet is working properly. While working, avoid introducing all the time objects inside the cabinet. This could prevent turbulences to the front barrier and to the laminar flow which would create abnormal situations and hinder a correct working, thus jeopardizing the cabinet functionality.

NEVER obstruct the front grid of the work surface not to compromise the protective barrier efficiency. Introduce in the working seat, the minimum quantity of necessary material, trying not to turn the cabinet in a stocking seat.

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THE WORKING SURFACE

The microbiological safety cabinets can be fitted with a working surface, perforated or partially closed.

In general, the perforated surface is to be preferred, for it is able to keep the laminar flow as most regular as possible, thus reducing the risk of cross-contaminations and obtaining a better product protection. In case instead of antitumoral drugs preparation, the closed surface turns out to be more comfortable for handling the phials and when using absorbent clothes. Though implying a better and more accurate manual skill to avoid the cross contaminations, this surface can be used all the way, given that the cross-contamination is not so important to each specific application.

MANUAL SKILL

A good manual skill is essential to properly work with a microbiological safety cabinet; the safety efficiency of the cabinet towards the product can be in fact reduced by the so-called "upstream contamination".

This phenomenon occurs each time the operator put a contaminated object (for example a hand) between the source of clean air (the HEPA filter) and the sterile manipulated material: the contaminants are transported by air through the non-sterile material and then conducted by the laminar flow itself towards the sterile object. The right way to operate wants that a sterile object is always as close as possible to the source of the laminar flow and that the material is moved through a parallel movement to the working surface and as a consequence perpendicular to the laminar flow. Only practice allows to acquire over time a proper manual skill.

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CLEANING AND MAINTENANCE OF THE MICROBIOLOGICAL SAFETY CABINET

A good manual skill allows to minimize the risk of drippings. In case of product leaks, it is necessary to promptly act in order to avoid contamination (sterilize and take off the working surface, sterilize the pouring through a proper bactericidal agent, eliminate the poured fluid by using the right absorbing media).

Day by day, after each work session, cleaning and disinfection of the cabinets and of the surface are required; use preferably some 70% ethyl alcohol or a disinfectant turned out to be efficient towards the used micro-organisms. Do not use chlorine-based compounds or disinfectants with chlorine because they could provoke phenomena such as rust on the stainless steel coatings of the cabinet.

Following the right daily cleaning procedure allows not to use the UV-rays germicide lamp. If you want to use the UV-radiation all the way, a thirty-minute radiation is sufficient to obtain the germicide effect. It is reminded that using the germicide lamp cannot replace at all the correct daily maintenance. It is expedient to replace the lamp every year or after 2000 hours of use.

Prepare some SOP (Standard Operating Procedures) as well as a programme of periodical maintenance aiming at monitoring and checking over time that the technical specifications declared by the manufacturer are observed. Twice a year (according to the EN ISO 14644-2) check the cleaning class and the laminar and barrier flow speed.

The life of the HEPA filters, by following the correct use conditions of the cabinet, is longer than 5000 working hours. At the end of this period, the filters must be replaced, preferably by qualified technicians. Before replacing them, it is necessary to sterilize the whole cabinet and the HEPA filters, by complying with what suggested by the manufacturer in the Instruction Manual.

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WASTE TREATMENT AND DISPOSAL

The HEPA filters and the biological material handled inside the microbiological safety cabinet are to be considered as hospital special and dangerous waste, that is why they have to be eliminated by abiding by the laws in force in this field.

USE OF CYTOSTATIC DRUG SAFETY CABINETS WHILE PREPARING THERAPIES WITH ANTINEOPLASTIC DRUGS

As for this special application, the WHO suggests to prefer Class II Biohazard vertical laminar flow cabinets, because they are able to contain the dispersion of drugs both in aerosol and in powder dispersed in the air. It is moreover expedient that this equipment is exclusively destined to the preparation of antitumoral therapies.

It is absolutely better not to use horizontal laminar flow cabinets, because the flow of contaminated air is directed towards the operator as well it is better not to use vertical laminar flow cabinets not belonging to the microbiological safety typology, for they are designed just to protect the handled product. Besides the international, European and national standards, above mentioned for the case of microbiological safety cabinets, it is wise that the cabinet complies with the German rule DIN 12980:2005 Zytostatika - Werkebanke, Anforderungen, Prufung, the only one at the moment, issued in Europe for the cytostatic drugs cabinet.

The requirements that cabinets should have to prove that this product is compliant with the rules must be the following: vertical laminar flow Class II microbiological safety cabinet provided with three HEPA filters. The main filter placed below the working surface, that comes into contact with 100% of potentially contaminated air must be easy to replace by a method which guarantees a proper protection for the operator and for the room where the chemiotherapies are arranged.

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
The cabinet should have a channel communicating outwards, if this is not easy to be carried out, we advise to put an activated charcoal filter on the expulsion seat.

Nowadays, it is still not clear whether the UV lamp is useful or not. According to some sources, it plays quite an important role in the degradation of some photosensitive drugs.

**Working Instructions
of a cytotoxic drug
Safety Cabinet**

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REFERENCES

- Federal Standard 209/e (1992)
- ISO 14644: 1 (1999); 2 (2000); 3 (2004); 4 (2001); 5 (2004); 6 (2004); 7 (2004)
- World Health Organisation: Laboratory Biosafety Manual - WHO (1983)
- U.S. Department of Health and Human Services, Center for Disease Control: Classification of etiologic agents on the basis of hazard - Atlanta (1976)
- U.S. Department of Health and Human Services, Center for Disease Control: Biosafety in microbiological and biomedical laboratories - Washington (1984)
- National Institutes of Health: NIH guidelines for recombinant DNA research - Bethesda (1979)
- British Standard Institution - BS 5726: Specifications for microbiological safety cabinets - London (1992)
- National Sanitation Foundation - NSF/ANSI 49: Standard for Class II (laminar flow) biohazard cabinet - Ann Arbor, Mich. (2002)
- Department of Health and Social Security: Code of Practice for the Prevention of Infection in Clinical Laboratories and Post-Mortem rooms - London (1978)
- Directive No. 2000/54/CE of September 18th 2000
- EN 12469: Biotechnology - Performance criteria for microbiological safety cabinets (2000)
- DIN 12980: Zytostatika - Werkbanke Anforderungen, Prufung (2005)
- Rules for a good preparation of drugs in pharmacies - NBP FU XI Ed.
- Categorisation of biological agents according to hazard and categories of containment - 4th Edition della ACDP (Advisory Committee on Dangerous Pathogens) 1995
- AS 2252: Biological Safety Cabinets (2004)
- AS 2567: Laminar flow cytotoxic drug safety cabinets
- AS 1807: Cleanrooms, workstations, safety cabinets and pharmaceutical isolators

BH-EN 2000

The **BH-EN 2000** microbiological safety cabinets have been designed to meet the requirements of the European Standard EN 12469:2000.

They are vertical laminar flow Class II biological safety cabinets, according to the definition reported by all world rules, with 70% of filtered air on an absolute H14 HEPA filter and then recirculated in laminar/unidirectional flow back in the working chamber and the remaining 30% exhausted after filtration on a H14 HEPA filter.

They are suitable for micro-organisms and pathogens manipulation, according to the international standards definitions, they are able to protect at the same time operator, product and environment. They can be used while handling pathogens during microbiology, virology, haematology, cell cultures, recombinant DNA applications, and so on. They have been

devised with a biodynamic tightness sloped plenum, thus guaranteeing all the pathogen material is kept inside the cabinet and preventing risks of filter bypass or contamination leakage within the laboratory.

Speed can be set automatically and run by a microprocessor to guarantee the right flow of laminar and protective barrier air, even in case of gradual HEPA filter obstruction.

They are available in BH-EN 2000 S version equipped with a motor fan or in BH-EN 2000 D version with 2 motor fans. Both the models are fitted with a managerial system with a microprocessor for monitoring all the functional parameters (laminar flow speed and front barrier/ environment and working room temperature/residual working hours of the HEPA filters, activated charcoal and UV lamp/ total working hours/ electrical black out signal/ HEPA filter obstruction level) with alphanumeric display. An audible /visual alarm system informs the user about any working anomaly.

The available sizes for the working surface are equivalent to about 90 cm (BH-EN 2003 S or D), about 120 cm (BH-EN 2004 S or D), about 150 cm (BH-EN 2005 S or D) and about 180 cm (BH-EN 2006 S or D).

The cabinet is provided with front laminated safety glass, which can be opened as an up-and-over door. On request, you could have optional fittings including support stand, drawer unit, UV lamp with safety closure, electronic Bunsen, etc.,.

The cabinet model BH-EN 2000 is certified and meets the requirements of the European Standards concerning the microbiological safety cabinets provided by LNE of Paris France for the NF EN 12469:2000 Standard and by TÜV Nord Cert of Hamburg, Germany for the DIN EN 12469:2000 Standard.

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The **BHG 2000** microbiological safety cabinets have been designed to meet the requirements of the European Standard EN 12469:2000.

They are vertical laminar flow Class II biological safety cabinets, according to the definition reported by all the world rules, with 70% of filtered air on an absolute H14 HEPA filter and then recirculated in laminar/unidirectional flow back in the working chamber and the remaining 30% expelled after filtration on a H14 HEPA filter.

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They are available in BHG 2000 S version equipped with a motor fan or in BHG 2000 D version with 2 motor fans. Both the models are fitted with a managerial system with a microprocessor for monitoring all the functional parameters (laminar flow speed and front barrier/environment and working room temperature/residual working hours of the HEPA filters, activated charcoal and UV lamp/total working hours/electrical black out signal/HEPA filter obstruction level) with alphanumeric display. An audible/visual alarm system informs the user about any working anomaly.

The available sizes for the working surface are of about 90 cm (BHG 2003 S or D), about 120 cm (BHG 2004 S or D) and about 180 cm (BHG 2006 S or D).

The cabinet is provided with electrically operated, vertically sliding safety glass sash. On request, you could have optional fittings including support stand, drawer unit, UV lamp with safety closure, electronic Bunsen, etc.,

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CYTOSAFE-N 2000

The **CYTOSAFE-N 2000** biological and cytotoxic safety cabinets have been designed to meet the requirements of the European Standard EN 12469:2000 concerning the microbiological safety cabinets and of the product standard DIN 12980:2005 concerning the cabinets destined to handle cytostatic drugs.

They are vertical laminar flow Class II biological safety cabinets, provided with three HEPA filters, according to the definition reported by all the world rules, with 70% of filtered air on an absolute H14 HEPA filter and then recirculated in laminar/unidirectional flow back onto the working chamber and the remaining 30% expelled after filtration on a H14 HEPA filter.

100% of the recirculated air is filtered again on an absolute H14 HEPA filter, set immediately below the working surface.

They are suitable for micro-organisms and pathogens manipulation, according to the international standards definitions and to the rules for cytostatic substances, they are able to protect at the same time operator, product and environment. They can be used while handling pathogens during microbiology, virology, haematology, cell cultures, recombinant DNA applications, but also for preparing antitumoral chemotherapies with cytotoxic and cytostatic drugs. They have been devised with a biodynamic tightness sloped plenum, thus guaranteeing all the pathogen and/or contaminated material is kept inside the cabinet and preventing risks of filter bypass or contamination leakage within the laboratory.

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The HEPA filter, placed below the working surface, can be replaced in an absolutely safe way for operators and environment without cabinet sterilisation when cytostatic substances are used inside.

They are fitted with a managerial system with a microprocessor for monitoring all the functional parameters (laminar flow speed and front barrier/environment and working room temperature/residual working hours of the HEPA filters, activated charcoal and UV lamp/total working hours/ electrical black out signal/HEPA filter obstruction level) with alphanumeric display. An audible/visual alarm system informs the user about any working anomaly.

The available sizes for the working surface are of about 90 cm (CYTOSAFE-N 2003), about 120 cm (CYTOSAFE-N 2004) and about 180 cm (CYTOSAFE-N 2006).

The cabinet model CYTOSAFE-N 2000 meets the requirements of the DIN 12980:2005 Standard (cabinets for cytostatic drugs).

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The **TWO 30** vertical laminar flow cabinets offer a Class 100 clean working area to work with cell cultures and for the applications in the microbiological field where non-pathogen biological material is used. They are moreover suitable to handle nucleic acids/thermo-cyclisations/DNA amplification. They are small-sized, about 80 cm, light and can run both in their overpressure version and in the negative pressure one, to the operator's discretion. Fitted with an HEPA absolute filter, with an efficiency better than 99.995% MPPS (H14) and with a synthetic-fibre pre-filter with a minimum efficiency of 80%-90% ASHRAE. Their lighting system is provided with a fluorescent tube.

On request, you could have optional fittings including support stand, drawer unit, UV lamp with safety closure, electronic Bunsen, etc.,

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The **BIO M** cabinets are Class 100 (M3.5) or ISO 5 vertical laminar flow cabinets which guarantee an excellent product protection, by offering a decontaminated working area thanks to the HEPA absolute filter with an efficiency better than 99.995% MPPS (H14). These cabinets can be used in industries as well as in laboratories, above all when handling cell cultures or non-pathogenic biological material. They have an epoxy powder painted steel structure.

Fully AISI-304 stainless steel structure is available as option.

The working surface is in stainless steel AISI 304L. The available sizes for this surface are of about 120 cm (Mod. BIO 48 M), about 150 cm (Mod. BIO 60 M) and about 180 cm (Mod. BIO 72 M).

There is the possibility to install an exhaust HEPA filter or an exhaust activated carbon filter, support stand, drawer unit, UV lamp, automatic regulation, electronic Bunsen, etc.

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The **MINILAB** vertical laminar flow cabinets are small working stations, about 50 cm, suitable for handling products that must be protected against non-pathogenic airborne contamination. They can be easily moved and positioned in different laboratories according to different needs. The cabinets are provided with a HEPA filter with an efficiency higher than 99.995% MPPS (H14) as well as with a synthetic-fibre filter.

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The **KBN** cabinets are Class 100 (M3.5) or ISO 5 horizontal laminar flow cabinets and they are suitable to keep the product in safety and to prevent external contaminations; the external construction is made of epoxy powder painted steel.

Fully AISI-304 stainless steel cabinets are available as option.

These cabinets are fitted with HEPA filters having an efficiency of more than 99.995% MPPS (H14) and a pre-filter whose minimum efficiency is 80%-90% ASHRAE, they have moreover a stainless steel working surface AISI 304L and a lighting system with fluorescent tubes. Available in four different versions, with a width of 90 cm, 120 cm, 150 cm and 180 cm.

On request, you could have optional fittings including support stand, drawer unit, UV lamp, automatic regulation, electronic Bunsen, closing curtains, etc.

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